

MIDLAND METROLOGY LIMITED

Mounting Sample Support Clip



Features:

- (1) For supporting samples in metallographic analysis process.
- (2) Ensure the sample in the mould has been fixed by the clip then add curing agent and epoxy resin or acrylic powder in proportion.
- (3) Make sure the sample is placed vertically. Convenient for in-depth scientific analysis of samples.
- (4) Three Types of Clip: Transparent Rolling Tyoe, SUS Rolling Type, Black Triangle Type.
- (5) Position the sample vertically to observe the cross section.

CLIP-2826P



Technical Specification

Color: Transparent
Material: Plastic
Package: 100PCS/Case

CLIP-28260



Technical Specification

Color: Transparent, Black, White
Shape: Ω
Material: Plastic
Package: 100PCS/Case

CLIP-2826S



Technical Specification

Material: Stainless Steel
Size: 22mm
Package: 100PCS/Case